

Use of a Closed Vial Milling System in Investigations of Carbon Black Ink Formulations: its Advantages and Disadvantages



Stephen R Bysouth† and Muhammad Ali*

†Automaxion SARL, Varengeuebec, France *N.E.D University of Engineering & Technology, Karachi, Pakistan



NED UNIVERSITY JOURNAL OF RESEARCH - APPLIED SCIENCES (2014) XI, 4, 15

Introduction to the problem

Carbon Black in Ink

What is carbon black?

Carbon black is a material produced by the incomplete combustion of heavy petroleum products and is a form of paracrystalline carbon that has a high surface-area-to-volume ratio.

Used as a printing ink pigment (surface area 30 - 138 m²/g)
Higher BET surface area carbon black also exist

How is the ink made conventionally?

Mixture of carbon black, dispersants, water.


These ingredients are dispersed together under impact and shear in a mill, usually by means of energetic beads (bead mill).

Ink chemistry by Chemistry World 1 March 2003
<https://www.chemistryworld.com/news/ink-chemistry/3002158.article>

A Conventional laboratory process

Lab scale re-circulating bead mill (here a Eiger Torrance Mini motormill 50).

- Sample volume
 - 70 mL (minimum)
- Conditions
 - 3750 RPM
 - 3 hours (determined by earlier experiments)




The problem with conventional lab techniques

Number of Ingredients x Number of Levels (concentrations) x nConditions (time, temperature, etc) = **Many experiments**

Many experiments:

- Consumes a lot of ingredients
- Produces a lot of waste
- Takes a long time
- A lot of cleaning




A milling system was sought that could:

Process many samples at the same time
Use smaller samples, but sufficient for testing (~10 ml)
Suitably disperse the sample!

Automaxion Mills

Based on Planetary Bead Milling



Milling jars replaced with vial holder

Mill for 10 ml in 'EPA' vials samples selected


Mills using 1, 2 or 4 adapters available



Experimental

Automaxion single position P6 mill selected

- 8 samples of 10 ml (in 20 or 30 ml vials)
- Adapters absorb/dissipate heat
- Glass vials mean no cleaning
- Samples can be stored in vials



Dispersing additives


High molecular weight non-ionic copolymeric dispersants

Dispersant	Supplier
• Solsperse 44000	• Lubrizol Corp
• BYK-190	• BYK-Chemie
• Tego 760W	• Evonik Tego

These were selected as candidate materials by quick 'add and settle' screening experiments

Screening experiments

Sinking test to determine the affinity of dispersant towards pigment



Pigment added to dispersant solution and left for 24 hours without agitation.

The opacity and volume of sediment were qualitatively analysed

Procedure to prepare pigment concentrates

- Step 1: Determine the theoretical amount of dispersant
- Step 2: Determine the maximum pigment loading
- Step 3: Optimise the dispersant amount
- Step 4: Analyse the dispersions

Example results for dispersion of Vulcan XC605 (mentioned as VXC605) carbon black pigment using Tego 760W (760W) as dispersant are presented here.

Analytical techniques

Particle size analysis. Indicates in combination:

- if agglomerates are broken up (milling)
- if sufficient dispersant is provided (titration)
- if dispersant is effective (chemistry)

Instruments

Malvern Zetasizer



TA instruments Rheometer



Rheology: A critical parameter for printing. Indicates:

- interactions between the particles, fluids and surfaces, e.g. charge, reversible building of structure, etc.

Storage stability (change in rheology over time). Indicates:

- Non-reversible build-up of structure (re-agglomeration)
- Migration of ingredients




Results

Results (Particle size analysis)

Dispersion name	Avg. particle size (nm)	
	Eiger mill	Automaxion mill
VXC605 22 wt%, 760W 12% DOWP	113.8	90.47
VXC605 22 wt%, 760W 15% DOWP	119.3	90.26
VXC605 22 wt%, 760W 18% DOWP	123.3	90.23

DOWP: Dispersant On the Weight of Pigment

Results (Particle size analysis)




Results (Zeta potential)


Dispersion name	Zeta potential (mV)	
	Eiger mill	Automaxion mill
VXC605 22 wt%, 760W 12% DOWP	-23.2	-17
VXC605 22 wt%, 760W 15% DOWP	-22.3	-29.5
VXC605 22 wt%, 760W 18% DOWP	-21	-33.4

Rheological characterisation


Percentage dispersant effect with different mills



Comparison of mills for two dispersant levels



Storage stability for the two mills



Preparation times (approx.)

Eiger: 15 minutes weighing ingredients and preparation

- 30 minutes pre-mixing
- 3 hours milling
- 30 minutes cleaning

Total = 4 hrs 15 mins per sample

3 samples = 12 hrs 45 minutes (or 34 hrs for 8 samples)

Automaxion: 15 minutes weighing ingredients and prep

- 0 minutes pre-mixing (not done)
- 4 hrs 30 mins milling for up to 8 samples at a time
- <2 minutes bead separation
- 0 minutes cleaning (not needed)

Total = 4 hrs 47 mins for only 1 sample

3 samples = 5 hrs 20 minutes (or 6 hrs 45 min for 8 samples)

Preparation times (reality)

Eiger

1 samples = 4 hrs 15 mins (or 34 hrs for 8 samples)

Reality = 1 sample per 8 hr day

Automaxion (P6)

8 samples = 6 hrs 45 min

Reality = 8 samples per 8 hr day

Or Automaxion (P5)

32 samples in day

Conclusions

- The mill is able to make pigment dispersions for printable inks
- Enabled rapid experiments
- Easy to use
- Produced representative inks
- Reproducible
- Some trends in particle size and zeta potential need further investigation
- Unable to process highly viscous samples but they were also difficult using a recirculating mill