Protein Stability and Aggregation Employing nanoDSF

The Prometheus Series

Dr. James Wilkinson James.wilkinson@nanotemper.co.uk

NanoTemper Technologies Ltd Sheraton House, Castle Hill, Cambridge, UK


NanoTemper Technologies GmbH Flößergasse 4, 81369 München, Germany

www.nanotemper-technologies.com



www.nanotemper-technologies.com







Sample format

24-capillary chips for automated handling 10 μl of sample / capillary

Same performance as Prometheus NT.48

- 20 points/minute per capillary (at 1 °C/min)
- Measure full sample set in 3 seconds
- Heat from 15°C to 95°C (110°C)
- Concentrations from 5 μg/ml to > 250 mg/ml
- Aggregation detection with backreflection
- New software for thermal and chemical unfolding

PR.ThermControl


PR.ChemControl

PR.TimeControl




Combinatory Evaluation of Thermal Unfolding and Aggregation






Simultaneous detection of conformational and colloidal stability

Decrease in pH destablizes antibody but prevents aggregation of the unfolded state



Thermal Stability of ADCs





nanoDSF Technology



Native conditions:

No fluorescent dyes, thus buffer independence, any detergent.

Ultra-high resolution:

Unmatched density of data points, more than 20 measurements per capillary and minute/°C.

Fast measurements:

Measure 24 / 48 capillaries "on-the-fly" within 3 seconds.

Broad application range:

Detect protein concentrations from $5 \mu g/ml$ to > 250 mg/ml.

Adress all your targets:

Soluble proteins, antibodies up to transmembrane or hydrophobic proteins.

Efficient research:

Easy to handle device, intuitive software with one-click-analysis and maintenance free instrument.



Prometheus

Come and discuss further at our exhibition booth and enter prize draw for free Franconian wine!