


# Malvern's Biopharmaceutical Solutions



Viscotek SEC – Separation-based **Orthogonal Analysis** of conformation and size of native protein and large aggregates

Archimedes RMM – **Subvisible particle quantification** and characterisation

Zetasizer range – Protein **size, charge, stability** and **rheology** analysis using DLS and ELS



# Viscosizer 200 – automated low volume viscosity screening

- › Microcapillary viscometry allows **automated high throughput viscosity** analysis
- › TDA allows accurate **sizing** in the nanometer and sub-nanometer range
- › **Viscosizer 200** - Rapid orthogonal formulation screening

